

## Introduction to OCTAVIUS 4D

**Dr. Jan U. Wuerfel**  
Project Manager, PTW-Freiburg

# Overview


1 Introduction

2 Nice features


3 The algorithm

4 Conclusion

# The OCTAVIUS Family




1999




2001



2002




2007



2010



2012



2012

# OCTAVIUS II vs. OCTAVIUS 4D



OCTAVIUS II



Fast setup



2D information



1 measurement plane only



OCTAVIUS 4D



4D information





No angular dependence





Setup needs longer

# OCTAVIUS II vs. OCTAVIUS 4D




# OCTAVIUS II vs. OCTAVIUS 4D

Which OCTAVIUS for IMRT, VMAT, RapidArc, TomoTherapy, Cyberknife, ...



# 3D $\gamma$ Analysis


## OCTAVIUS II and OCTAVIUS 4D



# 3D $\gamma$ Analysis




2D or 3D  $\gamma$  in a plane

OCTAVIUS II and OCTAVIUS 4D



Full Volume Analysis

OCTAVIUS 4D





# Overview

1 Introduction

2 Nice features

3 The algorithm

4 Conclusion



The screenshot displays the PTW-VenSoft software interface. The top-left panel shows a 2D DVH plot for 'DoseMap: ID: s\_HN1\_005' with a color scale from 0 to 150 mm. The top-right panel shows 'Results' statistics:


Statistics	
Number of Dose Points:	10,609
Evaluated Dose Points:	8,349 (78.7 %)
Passed:	8,330 (99.8 %)
Failed:	19 (0.2 %)
Result:	99.8 %

The bottom-left panel shows a 2D DVH plot for 'Measurement: Complete range'. The bottom-right panel shows a 3D RTView of a skull with a measurement plane. A red dashed box highlights the 'Compare RTView DVH' button at the bottom of the RTView panel.






- ▶ Calculated independent from TPS
- ▶ Calculation in patient anatomy
- ▶ Based on ion chamber measurements
- ▶ Fast calculation (approx. 2 min for VMAT plan)
- ▶ Optional software module
- ▶ Requires VeriSoft 6.0 and Rotation Unit




# Want to Change the Resolution?

*One system, choice of arrays*



## OCTAVIUS Detector 729


- ▶ 729 5 mm x 5 mm detectors
- ▶ 10 mm detector pitch
- ▶ 25 % field coverage
- ▶ Up to 27 cm x 27 cm field size



## OCTAVIUS Detector 1000 SRS

- ▶ 977 2 mm x 2 mm detectors
- ▶ 2.5 mm detector pitch
- ▶ Almost up to 100 % field coverage
- ▶ Up to 11 cm x 11 cm field size

# Wireless Inclinometer



- ▶ Bluetooth® 2.1
- ▶ Operating distance max. 10 m
- ▶ Operating time min. 16 h
- ▶ Charging time 7.5 h
- ▶ 3 rechargeable batteries (AA)
- ▶ Available: May 2014
- ▶ Upgrade kit (L981443)

# General Advantages

- ▶ Purely from measured data. No TPS input required
- ▶ The plan to be verified is not required for the verification
- ▶ No angular dependency because beam is always perpendicular
- ▶ Ionization chamber technology
  - ▶ No calibration at user site required
  - ▶ Long term stable
  - ▶ Reliable

# Overview

1 Introduction

2 Nice features

3 The algorithm

4 Conclusion

# The Algorithm (Abbreviated)

- ▶ At the current gantry angle: consider a detector of the array and the dose it has measured
- ▶ Construct a ray line through the current detector to the focus of the beam
- ▶ Determine the current field size from all irradiated detectors of the array and choose PDD accordingly
- ▶ Using the PDD and the dose of the chosen detector: reconstruct the dose along the ray line to the focus
- ▶ Do this for all detectors of the detector array
- ▶ Do this for all gantry angles
- ▶ Sum up all these dose values and sort them into 3D dose voxels



# The Algorithm (Full Version)

White Paper

October 2013  
D913.200.06/00  
**PTW**

## Dose reconstruction in the OCTAVIUS 4D phantom and in the patient without using dose information from the TPS

B. Allgaier, E. Schüle, J. Würfel

PTW-Freiburg Physikalisch-Technische Werkstätten Dr. Pöchlau GmbH  
Lörracher Straße 7, 79115 Freiburg, Germany

### 1. Introduction

Intensity modulated radiation therapy (IMRT) and volumetric modulated arc therapy (VMAT) are state-of-the-art irradiation techniques for the delivery of highly conformal radiation fields to the target volume. These techniques require complex treatment planning system (TPS) algorithms as well as sophisticated irradiation methods. As a result, the use of quality assurance tools for the verification of the planned dose distribution prior to the treatment of the patient has become a standard procedure in clinical routine.

Early quality assurance (QA) tools for this purpose are two-dimensional devices based on stationary 2D detector panels [1], [2], [3], [4], [5]. They allow the measurement of one plane of dose values in order to compare it with the corresponding plane calculated by the TPS. Some of these tools have restrictions regarding the angle of the incident beam because of the angular dependence of the detector panel [4], requiring so-called collapsed beam measurements during which the gantry of the treatment unit is fixed at a certain angle. Other tools allow composite plan measurements by correcting for the directional dependence either by the application of correction factors [1], [5], or by the use of suitably shaped phantoms [3].

Because of the complex volumetric shape of the conformal fields the limitation to one measuring plane is considered a disadvantage. Ideally, a full 3D dose matrix is desirable. To date the only realistic 3D dosimetry system is based on gel dosimetry [6], but unfortunately there are restrictions in the use of such systems in clinical routine such as the availability of gels with reproducible features, temperature dependence, the change of the irradiated gel volume by time as well as the necessity of relatively complex equipment for the evaluation of the irradiated gel samples.

measurements, combined with a dose reconstruction method that determines dose values in the 3D volume. Such 2D detectors are available either as planar detector panels [7], [8], [9] or as the surface of a cylinder [10]. All these systems measure time-resolved dose values at a limited number of positions and reconstruct 3D dose values in the complete volume. While some systems [7], [9] perform this task as independent measuring tools other systems [8], [10] require the dose reconstruction. Many clinical physicists prefer QA tools to be independent of the object to be tested. Therefore, the OCTAVIUS 4D was designed as a truly independent tool for pre-treatment quality assurance. This paper describes the algorithms on which the OCTAVIUS 4D dose reconstruction in the phantom and in the patient are based. In addition, some of the results of the system verification are presented.

### 2. 3D dose reconstruction methods

#### 2.1 TPS dependent methods

One example of a QA system for 3D dose measurements is the ArcCHECK with 3DVH option [10]. The ArcCHECK detectors are located in a plane that forms the surface of a cylinder, allowing the radiation beam to hit at least the central detectors perpendicularly at any gantry angle. The reconstruction of the dose matrix calculated by the TPS. Therefore, before ArcCHECK can determine the volumetric dose grid, the treatment plan including its dose values needs to be entered into the ArcCHECK software. Then the TPS data are "corrected" by the measured values, and the corrected TPS data are presented as the "ArcCHECK measuring result". This method is called "ArcCHECK Planned Dose Perturbation" (ACPPD) [10]. The scaling factors for the perturbation of the plan are interpolated from the measured entry and exit values of the phantom, resulting in a "morphing" of the dose grid rather than a simple re-normalization in order to determine the

# Overview

1 Introduction

2 Nice features

3 The algorithm

4 Conclusion

# Take Home Message OCTAVIUS 4D

- ▶ Purely measured data, no TPS dose input required
- ▶ DVH 4D: also purely measured data, no TPS dose input required
- ▶ The plan to be verified is not required for the verification
- ▶ Understandable algorithm, no black-box technology
- ▶ Choice of arrays available
- ▶ Reliable technology:
  - ▶ No angular dependence
  - ▶ No user-calibration necessary
  - ▶ Long term stable

# Thank you for your attention

